



CardioChek P-A™

Point of Care Test System

Professional

[Overview](#)
[Products](#)
[Technology](#)
[Training](#)
[Accuracy](#)
[In-Office Testing](#)

Test NOW

Know NOW


Act NOW

CardioChek P-A Technology

The underlying technology of the CardioChek P-A System employs innovative concepts to emulate the complex chemical analysis, sample preparation and measurement methods of clinical labs in a simple format suitable for untrained operators.

Operating Principle

A precise volume of whole blood is dispensed on to the application window of a PTS Panels Test Strip that has been inserted into a CardioChek P-A analyzer. Red cells and other non-plasma components are removed from the plasma. The plasma reacts in a precisely controlled manner with a series of chemically impregnated membranes, ultimately producing a color change in the final membrane layer. The analyzer measures the color intensity from the color spectral profile and compares to lot-specific calibration information contained in the MEMo Chip to calculate an accurate result. The analyzer then displays the result, simultaneously triggers a printout (if the dedicated printer is attached), and stores the test information in memory for later recall.



CardioChek P-A™ Analyzer

The CardioChek P-A is a state-of-the-art, hand-held, battery-operated, software-driven reflectance spectrophotometer.

PTS Panels Test Strips

PTS Panels Test Strips are composed of several layers of reagent-impregnated membranes sandwiched to a plastic carrier. Each membrane layer performs specific functions to process and facilitate chemical reactions. The phase chemistry membranes are typically stable at room temperatures for 18 to 24 months.

MEMo Chip

MEMo Chips are individualized for each lot of PTS Panels Test Strips and contain calibration settings, expiration dating and other information that calibrates the CardioChek P-A analyzer to every strip it "reads." A MEMo Chip is packaged with every vial of strips to ensure accurate calibration specific to each lot of strips.



Benefits of CardioChek P-A Technology

- **Accurate.** Factory calibrated with MEMo Chip
- **Fast.** Results in approximately 1-2 minutes*
- **Versatile.** Ability to test multiple-test panels or individual tests
- **Convenient.** Small sample (15-40µL, fingerstick or venous (EDTA or heparin))
- **Reliable.** No moving parts
- **Easy, simple.** Stick finger > collect sample > dispense sample > read results
- **CLIA-Waived.***
- **Stable.** Room temperature strip storage 18-24 months (no refrigeration required*)
- **Portable.** Test anywhere, any time (pocket-size, battery-operated)
- **Inexpensive.**
- **Profitable.** Covered by most US payors for medically-indicated testing
- **Memory.** Stores last 30 tests for recall
- **Printer and Computer Connectivity.** Optional

*Creatinine requires refrigeration, provides test results in 8 minutes, and is not CLIA-Waived.

Test NOW

Know NOW

Act NOW




Customer Service

[Overview](#)
[Professional Training](#)
[Consumer Training](#)
[Frequently Asked Questions \(FAQ\)](#)
[Help Hotline](#)
[Quality Control](#)
[Control Solution Range Values](#)
[Documentation](#)
[Specifications](#)
[Professional Catalog](#)
[Consumer Catalog](#)
[Medical Distributors](#)
[Consumer Distributors](#)

Product Specifications

CardioChek P•A
 CardioChek
 PTS Panels Test Strips
 CardioChek P•A Printer

CardioChek Analyzers

CardioChek P•A:

Length: 5.5 in. (14.0cm)
 Width: 3.0 in. (7.6 cm)
 Height: 1.0 in. (2.5 cm)
 Weight: 4.3 oz. (121.9 g)
 Requires 2 AAA batteries
 Battery life: approximately 300 tests
 CLIA: Waived

CardioChek:

Length: 5.5 in. (14.0cm)
 Width: 3.0 in. (7.6 cm)
 Height: 1.0 in. (2.5 cm)
 Weight: 4.3 oz. (121.9 g)
 Requires 2 AAA batteries
 Battery life: approximately 300 tests
 CLIA: Waived

PTS Panels Multiple Test Strips

Lipid Panel:

Range:	Total Cholesterol:	100-400 mg/dL (2.59-10.36 mmol/L)
	HDL Cholesterol:	15-100 mg/dL (0.39-2.59 mmol/L)
	Triglycerides:	50-500 mg/dL (0.56-5.65 mmol/L)
Time:		about 2 minutes
Sample:		35-40 μ L whole blood (fingerstick or venous)
Analyzer:		CardioChek P•A only
CLIA:		Waived

Cholesterol + Glucose:

Range:	Total Cholesterol:	100-400 mg/dL (2.59-10.36 mmol/L)
	Glucose:	20-600 mg/dL (1.11-33.3 mmol/L)
Time:		about 2 minutes
Sample:		25-30 μ L whole blood (fingerstick or venous)
Analyzer:		CardioChek P•A only
CLIA:		Waived

Cholesterol + HDL:

Range:	Total Cholesterol:	100-400 mg/dL (2.59-10.36 mmol/L)
	HDL Cholesterol:	15-100 mg/dL (0.39-2.59 mmol/L)
Time:		about 2 minutes
Sample:		35-40 μ L whole blood (fingerstick or venous)
Analyzer:		CardioChek P•A only

PTS Panels Individual Test Strips**Total Cholesterol:**

Range: 100-400 mg/dL (2.59-10.36 mmol/L)
Time: about 2 minutes
Sample: 15 μ L whole blood (fingerstick or venous)
Analyzer: CardioChek P•A or CardioChek
CLIA: Waived

HDL Cholesterol:

Range: 25-85 mg/dL (0.65-2.20 mmol/L)
Time: about 1 minute
Sample: 15 μ L whole blood (fingerstick or venous)
Analyzer: CardioChek P•A or CardioChek
CLIA: Waived

LDL Cholesterol:

Range: 50-200 mg/dL (1.3-5.18 mmol/L)
Time: about 2 minutes
Sample: 15 μ L whole blood (fingerstick or venous)
Analyzer: CardioChek P•A only
CLIA: Waived

Ketone:

Range: 2-70 mg/dL (0.19-6.72 mmol/L)
Time: about 2 minutes
Sample: 15 μ L whole blood (fingerstick or venous)
Analyzer: CardioChek P•A or CardioChek
CLIA: Waived

Triglycerides:

Range: 50-500 mg/dL (0.56-5.65 mmol/L)
Time: about 1 minute
Sample: 15 μ L whole blood (fingerstick or venous)
Analyzer: CardioChek P•A or CardioChek
CLIA: Waived

Glucose:

Range: 20-600 mg/dL (1.11-33.3 mmol/L)
Time: about 1 minute
Sample: 15 μ L whole blood (fingerstick or venous)
Analyzer: CardioChek P•A or CardioChek
CLIA: Waived

Creatinine:

Range: 0.2-10 mg/dL (17.8-885 μ mol/L)
Time: within about 8 minutes
Sample: 20 μ L whole blood (fingerstick or venous)
Analyzer: CardioChek P•A
CLIA: Moderate Complexity

CardioChek P•A Printer

Input: 100-240 VAC, 47-63 Hz, 1.5 A

Output: +20vdc, 2.5A

Length: 7.0 in. (17.78 cm)

Width: 3.5 in. (8.89 cm)

Height: 5.5 in. (13.97 cm)

Weight: 2.0 lbs. (0.97 kg)

Printer Label Specifications: Size: 2 ¼" x 3" per label, 500 labels per roll

Printer Paper Specifications: Size: 2 ¼" x 200 ft. roll

[SITE MAP](#) | [PRODUCTS](#) | [QUALITY CONTROL](#) | [PROFESSIONAL](#) | [CONSUMER](#) | [CONTACT PTS](#)

© 2006 Polymer Technology Systems, Inc.



CardioChek P-A™

Point of Care Test System

Professional

Overview

Products

Technology

Training

Accuracy

NCEP Guidelines

Proficiency Surveys

Technical Documentation

In-Office Testing

Test NOW

Know NOW

Act NOW

NCEP Guidelines

Total System Error

The generally accepted standard for quality of a result from any test method relates to the difference between the result obtained and the reference value. That difference is referred to as "Total System Error", which is a measure of both accuracy of the system (i.e., the bias of the system, or how close the average of a series of results on the same sample are to the reference value) and precision (i.e., the reproducibility of the system, or how close a series of results on the same sample are to each other).

Reference Value

Because there is variability and system error in all commercial methods for determining lipid concentration, determining the reference value can be difficult. The reference value is typically defined as the result given by the standard methods established by the Center for Disease Control and Prevention (CDC) for lipid analyses (the "National Reference System for Cholesterol" (NRS/Chol)).

CRMLN Reference Labs

The CDC established reference laboratories in the US ("Cholesterol Reference Method Laboratory Network" (CRMLN)) to evaluate and certify manufacturers' lipid products against reference standards.

NCEP Guidelines for Total System Error

The National Institutes of Health (NIH), through its National Cholesterol Education Program (NCEP), has established test protocols and guidelines for acceptable deviation from "truth", or the NRS/Chol reference value. These guidelines state that Total System Error (TSE) should be within the following limits from reference when these test protocols are followed:

Analyte	Total Error
Total Cholesterol	±8.9%
HDL Cholesterol	±13%
Triglycerides	±15%
LDL Cholesterol (direct measured)	±12%

Test NOW

Know NOW

Act NOW



Professional

Overview

Products

Technology

Training

Accuracy

NCEP Guidelines

Proficiency Surveys

Technical Documentation

In-Office Testing

Test NOW

Know NOW

Act NOW

CardioChek P•A™
Point of Care Test System

Proficiency Surveys

Proficiency Survey Tests

CLIA regulations require that all blind moderately-complex and highly-complex clinical laboratories (and in some even waived laboratories) regularly participate in "Proficiency Surveys" using unknown specimens to monitor accuracy and precision of results obtained by the methods and products they utilize. There are several US and international proficiency testing programs, including surveys conducted by the American Proficiency Institute (API) and the College of American Pathologists (CAP).

Standards for Acceptable Results

In these proficiency surveys the target value is usually regarded as the average of results submitted by all participants of the survey across all methods and products. Lipid results that are within ± 2 standard deviations of the mean analyzer type used are usually considered acceptable.

Examples of Accuracy and Precision of Common Commercial Lipid Test Methods

In the 2004, 3rd Test Event proficiency survey conducted by the American Proficiency Institute (See API's Part 1 Data Summaries) with 2,636 participants, even large automated analyzers reported significant bias from the mean value and significant variation in results, even within labs using the same methodologies. For example:

<u>Total Cholesterol</u>	<u>Mean</u>	<u>Bias*</u>	<u>SD</u>	<u>CV*</u>	<u>Range</u>
All Participants	195.1		8.7	4.5%	175-215
Beckman Synchron CX	194.2	+0.5%	5.3	2.7%	174-214
Cholestech LDX	203.8	+4.5%	7.4	3.6%	189-219
Dade Dimension	186.5	-4.4%	4.3	2.3%	167-206
OCD Vitros	199.6	+2.3%	6.1	3.1%	179-220
Roche Hitachi	198.3	+1.6%	7.4	3.7%	178-219

Bias is defined herein as the difference between the mean result for the test system and the mean result for all participants.
CV is defined herein as the Standard Deviation as a percent of the Mean.

<u>HDL Cholesterol</u>	<u>Mean</u>	<u>Bias*</u>	<u>SD</u>	<u>CV*</u>	<u>Range</u>
All Participants	45.3		8.1	17.9%	31-59
Beckman Synchron CX	54.8	21.0%	4.0	7.3%	38-72
Cholestech LDX	43.6	-3.8%	7.8	17.9%	28-60
Dade Dimension	46.7	3.1%	6.2	13.3%	32-61
OCD Vitros 250-950	39.7	-12.4%	3.8	9.6%	27-52
Roche Hitachi	34.9	-23.0%	2.2	6.3%	24-46

<u>Triglycerides</u>	<u>Mean</u>	<u>Bias*</u>	<u>SD</u>	<u>CV*</u>	<u>Range</u>
All Participants	150.6		9.0	6.0%	112-189
Beckman Synchron CX	150.5	-0.1%	8.2	5.1%	112-189
Cholestech LDX	159.5	5.9%	8.2	5.1%	143-176
Dade Dimension	153.6	2.0%	5.8	3.8%	115-192
OCD Vitros 250-950	159.0	5.6	6.0	3.8%	119-199
Roche Hitachi	154.8	2.8%	3.9	2.5%	116-194

<u>Direct LDL</u>	<u>Mean</u>	<u>Bias*</u>	<u>SD</u>	<u>CV*</u>	<u>Range</u>
All Participants	87.1		20.2	23.2%	46-128
Beckman Synchron CX	64.4	-26.1%	4.2	6.5%	56-73
Cholestech LDX	Not Available				
Dade Dimension	84.3	-3.2%	5.4	6.4%	73-96
Roche Hitachi	111.1	27.6%	2.9	2.6%	105-117

CardioChek P•A System Proficiency

Launched in mid-2003 the CardioChek P•A System was too new to participate in 2003 and 2004 proficiency surveys.

will, however, be listed in 2005. In the meantime, data submitted to the FDA as part of the 510(k) clearance procedure demonstrate that the CardioChek P•A with Lipid Panel strips and the new direct LDL strips perform comparably above survey statistics:

	<u>Mean</u>	<u>SD</u>	<u>CV</u>
Total Cholesterol*	197.2	8.4	4.3%
HDL-Cholesterol*	39.2	2.5	6.4%
Triglycerides*	157.0	6.1	3.9%
LDL-Cholesterol, Direct**	113.6	6.11	5.4%

* Data from Lipid Panel package insert.

** Data from LDL package insert.

TestNOW

KnowNOW

ActNOW

[SITE MAP](#) | [PRODUCTS](#) | [QUALITY CONTROL](#) | [PROFESSIONAL](#) | [CONSUMER](#) | [CONTACT PTS](#)

© 2006 Polymer Technology Systems, Inc.